

Fine mist. Foldable form.

Foldable eye-misting eyewear

An eyewear format for mist delivery near the eyes, paired with the specified solution.



Wearing-scene photograph

Wearable format

Leaves the hands free during stationary use.

Foldable temples

Folds for convenient storage between uses.

Removable parts

Detachable reservoirs and removable lens assemblies.

Measured operation: about 1.5-2 hours


Per charge; operating time depends on the selected mode.

Battery runtime is not a recommended continuous wearing duration.

For optical retailers and distribution partners.

See the mist. Explore the design.

Actual device operation and removable components support a clear product review.



Actual operation of the supplied unit



Reservoir, lens and frame assemblies

Mist delivery	Modules sit beside the lenses. The companion video shows actual mist operation.
Removable parts	Reservoir and lens assemblies can be removed. Confirm the supplied configuration.
Lens options	Discuss prescription-lens configuration and fitting requirements with your partner.
Specified solution	For use with the specified sodium-hyaluronate misting solution.

From a sample review to a channel fit

Assess fit, mist sensation, visibility and component handling, then discuss configuration and commercial terms.

Individual wearing experience varies. For product-selection discussion; local classification and sales requirements need review before launch.